

상세규격서

구분	품명	단위	수량	비고
1	1. Right-Angle Kinematic Mirror Mount with Tapped Cage Rod Holes 60 mm Cage System and SM2 Compatible M4 and M6 Mounting Holes(Model: KCB2/M)	ea	10	
2	2. Right-Angle Kinematic Mirror Mount with Smooth Cage Rod Bores 60 mm Cage System and SM2 Compatible M4 and M6 Mounting Holes (Model: KCB2C/M)	ea	10	
3	3. 60 mm Cage Plate SM2 Threads 0.5" Thick M4 Tap (Two SM2RR Retaining Rings Included) (Model: LCP34/M)	ea	20	
4	4. 60 mm Cage Plate SM2 Threads 0.9" Thick M4 Tap (Two SM2RR Retaining Rings Included) (Model: LCP34T/M)	ea	10	
5	5. 60 mm Cage Alignment Plate with Ø1 mm Hole (Model: LCPA1)	ea	5	
6	6. Snap-On 60 mm Cage Mounting Bracket M4 Tap (Model: LCPMA4/M)	ea	10	
7	7. Vertical 30 mm Cage Clamp for Ø1.5" Posts Metric (Model: CV1530/M)	ea	8	
8	8. 30 mm Cage Clamp for Ø1.5" Posts Included Quick-Release Handle Metric (Model: CH1530/M)	ea	5	
9	9. 60 mm Cage Clamp for Ø1.5" Posts Included Quick-Release Handle Metric (Model: CH1530/M)	ea	5	
10	10. M4 x 0.7 Alloy Steel Nylon-Tipped Setscrew 4 mm Long 10 Pack (Model: SS4MN4)	ea	10	
11	11. Magnetic Quick-Release Carriage Set M4 Tap (Model: KB1P/M)	ea	20	
12	12. Kinematic Prism Mount 25.4 mm Deep M4 Taps (Model: KM100PM/M)	ea	10	
13	13. Kinematic Prism Mount 54.4 mm Deep M4 Taps (Model: KM200PM/M)	ea	10	

상세내역

1. 구성

1. Right-Angle Kinematic Mirror Mount with Tapped Cage Rod Holes 60 mm Cage System and SM2 Compatible M4 and M6 Mounting Holes(Model : KCB2/M)
2. Right-Angle Kinematic Mirror Mount with Smooth Cage Rod Bores 60 mm Cage System and SM2 Compatible M4 and M6 Mounting Holes (Model : KCB2C/M)
3. 60 mm Cage Plate SM2 Threads 0.5" Thick M4 Tap (Two SM2RR Retaining Rings Included) (Model : LCP34/M)
4. 60 mm Cage Plate SM2 Threads 0.9" Thick M4 Tap (Two SM2RR Retaining Rings Included) (Model : LCP34T/M)
5. 60 mm Cage Alignment Plate with Ø1 mm Hole (Model : LCPA1)
6. Snap-On 60 mm Cage Mounting Bracket M4 Tap (Model : LCPMA4/M)
7. Vertical 30 mm Cage Clamp for Ø1.5" Posts Metric (Model : CV1530/M)
8. 30 mm Cage Clamp for Ø1.5" Posts Included Quick-Release Handle Metric (Model : CH1530/M)
9. 60 mm Cage Clamp for Ø1.5" Posts Included Quick-Release Handle Metric

(Model : CH1560/M)

10. M4 x 0.7 Alloy Steel Nylon-Tipped Setscrew 4 mm Long 10 Pack (Model : SS4MN4)
11. Magnetic Quick-Release Carriage Set M4 Tap (Model : KB1P/M)
12. Kinematic Prism Mount 25.4 mm Deep M4 Taps (Model : KM100PM/M)
13. Kinematic Prism Mount 54.4 mm Deep M4 Taps (Model : KM200PM/M)

2. 구성

1. Right-Angle Kinematic Mirror Mount with Tapped Cage Rod Holes 60 mm Cage System and SM2 Compatible M4 and M6 Mounting Holes (Model : KCB2/M)



- Right-Angle Kinematic Mirror Mount
 - Tapped Cage Rod Holes
 - 60 mm Cage System and SM2 Compatible
 - M4 and M6 Mounting Holes
2. Right-Angle Kinematic Mirror Mount with Smooth Cage Rod Bores 60 mm Cage System and SM2 Compatible M4 and M6 Mounting Holes (Model : KCB2C/M)



Back




Front

KCB2C

- Right-Angle Kinematic Mirror Mount
- Smooth Cage Rod Bores
- 60 mm Cage System and SM2 Compatible
- M4 and M6 Mounting Holes

3. 60 mm Cage Plate SM2 Threads 0.5" Thick M4 Tap (Two SM2RR Retaining Rings Included) (Model : LCP34/M)

LCP34



Accepts Ø2" Optics

- 60 mm Cage Plate
- SM2 Threads
- 0.5" Thick
- M4 Tap (Two SM2RR Retaining Rings Included)

4. 60 mm Cage Plate SM2 Threads 0.9" Thick M4 Tap (Two SM2RR Retaining Rings Included) (Model : LCP34T/M)

LCP34T



- 60 mm Cage Plate
- SM2 Threads
- 0.9" Thick
- M4 Tap (Two SM2RR Retaining Rings Included)

5. 60 mm Cage Alignment Plate with $\text{\O}1$ mm Hole (Model : LCPA1)



- 60 mm Cage Alignment Plate
- $\text{\O}1$ mm Hole

6. Snap-On 60 mm Cage Mounting Bracket M4 Tap (Model : LCPMA4/M)



- 60 mm Cage System Compatible
- Snap-On Functionality for Easy Installation or Removal
- Firmly Lock onto $\text{\O}6.0$ mm ER Rods with Flexure Slot Locking Design
- Off-Axis Mounting Adjustment Range
- LCPMA3: ± 0.88 " (± 22.3 mm)
- LCPMA4(M): ± 0.75 " (± 19.1 mm)
- Brackets Offer Clearance for SM2 Lens Tubes

7. Vertical 30 mm Cage Clamp for $\text{\O}1.5$ " Posts Metric (Model : CV1530/M)



CV1530

- Mount 30 mm Cage Systems to Our $\text{\O}1.5$ " Posts
- Can Be Oriented to Mount Cage Systems Vertically or Horizontally
- Quick-Connect Mechanism for Fast Drop-In/Drop-Out of Cage Assemblies
- 15 lb(6.8 kg) Load Capacity

8. 30 mm Cage Clamp for $\text{\O}1.5$ " Posts Included Quick-Release Handle Metric (Model : CH1530/M)



- Securely Mount 30 mm Cage Systems to our Ø1.5" Posts
- Quick-Connect Mechanism for Fast Drop-In/Drop-Out of Cage Assemblies using Ø6 mm Cage Rods
- Cage Clamping Mechanism Activated by 3/16" (5 mm) Hex Cap Screw
- Includes Quick-Release Handle for Easy Mounting to Ø1.5" Posts

9. 60 mm Cage Clamp for Ø1.5" Posts Included Quick-Release Handle Metric (Model : CH1560/M)



- Securely Mount 60 mm Cage Systems to our Ø1.5" Posts
- Quick-Connect Mechanism for Fast Drop-In/Drop-Out of Cage Assemblies using Ø6 mm Cage Rods
- Cage Clamping Mechanism Activated by 3/16" (5 mm) Hex Cap Screw
- Includes Quick-Release Handle for Easy Mounting to Ø1.5" Posts

10. M4 x 0.7 Alloy Steel Nylon-Tipped Setscrew 4 mm Long 10 Pack (Model : SS4MN4)



- M4 x 0.7 Alloy Steel Nylon-Tipped Setscrew
- 4 mm Long
- 10 Pack

11. Magnetic Quick-Release Carriage Set M4 Tap (Model : KB1P/M)



- SM1-Threaded (1.035"-40) Mounting Carriage for Ø1" Optics
- One Retaining Ring (SM1RR) Included
- Magnetically Coupled Removable Mounting Carriage
- 8-32 (M4) Tapped Hole on Back Plate for Post Mounting
- Ideal for Use with Lenses, Filters, Apertures, and Other Optics

12. Kinematic Prism Mount 25.4 mm Deep M4 Taps (Model : KM100PM/M)



- Two 80 or 100 TPI Adjusters for $\pm 4^\circ$ of Fine Tip and Tilt Control
- Platform can be Moved to Either Arm for Left- or Right-Handed Configuration
- Designed for Use with Clamping Arms (Sold Below)
- #8 (M4) Through Holes or Taps in Back Plate of Mounts for $\varnothing 1/2$ " Post Mounting
- Fabricated from Black Anodized Aluminum

13. Kinematic Prism Mount 54.4 mm Deep M4 Taps (Model : KM200PM/M)



- Two 80 or 100 TPI Adjusters for $\pm 4^\circ$ of Fine Tip and Tilt Control
- Platform can be Moved to Either Arm for Left- or Right-Handed Configuration

- Designed for Use with Clamping Arms (Sold Below)
- #8 (M4) Through Holes or Taps in Back Plate of Mounts for Ø1/2“ Post Mounting
- Fabricated from Black Anodized Aluminum

3. 기타사항

가. 설치위치 : 경북 구미시 구미대로 350-27 XR 디바이스 개발 지원센터 104호

나. 설치개요 : 규격서의 조건에 맞는 광학 물품 납품

다. 납품기한 : 계약일로 부터 60일 이내

라. 계약의 이행

- 1) 본 규격서에 대한 주계약자(이하 “계약자” 라 한다)와 발주기관 간의 해석상 이견은 발주기관의 해석에 따라야 하며, 본 규격서에 특별히 명시되지 않은 사항은 정부계약법 등의 규정과 일반 관례에 따른다.
- 2) 발주기관은 계약내용 전부에 대한 권리를 행사할 수 있으며 계약자는 이에 대하여 전적으로 동의해야 한다.
- 3) 본 규격서에 명시된 모든 조항은 본 사업의 원활한 추진에 필요한 최소한의 사항만을 규정하였으므로 상세히 기술하지 않았거나 누락된 사항에 대해서도 정상운영 및 서비스에 문제가 발생되지 않도록 계약자는 사전 조치를 취하여야 한다.
- 4) 본 규격서에서 요구하는 성능 및 특성에 대해 계약자의 오인 또는 누락된 상태로 계약 체결되어, 계약 후 이를 보완하기 위한 설계 변경, 물품 등을 추가해야 할 경우 계약자 부담으로 시행한다.
- 5) 계약자는 공급하는 제품의 성능, 품질 등에 대한 최종적인 책임을 진다.

4. 기타조건

가. 납품기한 : 계약 후 60일 이내

다. 규격관련 문의사항 : 강필성 선임 : 010-4261-0006